


AMRAD  **Medical**



AAU *Classic*


U-ARM RADIOGRAPHIC SYSTEM



Large display of arm and detector rotation



Large LCD display with easy touch buttons for accurate positioning




Near-patient controls offer efficient equipment positioning refinements



Position equipment from anywhere in the room



Choose from a variety of mobile tables to add to your system. Available, if needed.



Peak Performance

Achieve peak performance by migrating most exams to standing positions improving productivity and throughput.





Smart Investment

AAU Classic offers equal features and earn equal reimbursement to all of its more costly competition and all of this at 5 times the warranty. Smart decisions equal smart investments. Be smart!



Refine equipment positions and remove complicated equipment angles with conveniently located control panel

**Add the AMRAD
Medical DR Elite for
additional features
and automation.**



AMRAD  **Medical**

www.amradmedical.com

SUMMIT INDUSTRIES, LLC
 7555 N. Caldwell Avenue, Niles, Illinois 60714
 (p) 773-353-4000 | (p) 1-800-729-9729

Certified System

 SAI GLOBAL
 ISO 13485
 Quality



LITR040-1120