

DELWORKS EDR

E14Ce Wireless Detector

Affordable and Technology Rich

Features

Light, 6.8 lbs portable 14 x 17"

Internal accelerometer that senses motion, in and out of ready/sleep modes, extending battery life up to 8 hours

EasyConnect
Automatic Exposure Detection (AED)




E14Ce battery and charging dock

Specifications

Detector Size	14" x 17" (35 x 43 cm) ISO 4090 Cassette Size
Dimensions	15.1" x 18.1" x 0.6" (383.5 x 459.5 x 15 mm)
Active Area	13.8" x 16.9" (350 x 430 mm)
Type	Wireless
Scintillator	Cesium Iodide (CsI) <small>with Amorphous Silicon Photodiode (a-Si)</small>
Pixel Pitch	160 μ m
Limiting Resolution	3.5 lp/mm
Pixel Matrix	2400 x 2880
Bit Depth (A/D Conversion)	16-bit
Image Display	>6 seconds
Cycle Time	9 seconds
Exposure Window	Up to 3.2 seconds
MTF	61% @ 1 lp/mm 30% @ 2 lp/mm 14% @ 3 lp/mm
DQE	70% @ 0 lp/mm 51% @ 1 lp/mm 42% @ 2 lp/mm 22% @ 3 lp/mm
Weight of battery and detector	6.8 lbs. (3.1 kg)
Battery Type	Li-ion 3.35Ah 7.4V
WiFi	802.11a/b/g/n, 2.4GHz/5GHz
Drop Sensor(s)	Mechanical

15.1" (383.5 mm)

0.6" (15 mm)



18.1" (459.5 mm)

11.2018 #9000-DW-E14Ce



www.delmedical.com



New York
28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Chicago
241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

DEL MEDICAL
A **UMG** COMPANY