

DENTAL IMAGING EQUIPMENT YOU CAN TRUST





The versatility offered by four distinct wall mounting options allows unprecedented mounting flexibility.

Image-Vet™ DC

The Image-Vet DC is our premier x-ray generator featuring superior image quality and is seamlessly compatible with the EVA Vet Digital Dental Imaging System. Advanced DC technology results in a lighter tubehead and more stable positioning experience. The unique wireless remote allows you to program and take exposures from a distance.

➡ **FREE DC-101 AND POSITIONING ONLINE TRAINING AVAILABLE WITH EVERY IMAGE-VET DC PURCHASE**



Unique DCR wireless remote allows the operator to set time and exposure and take the image from a distance. Remote also features manual and anatomical programming.



Optional mobile stand gives you ultimate flexibility. Features 4" medical-grade casters and dual-locking brakes for optimal mobility and stability.